ES

# EIA CanAg NSE

REF

420-10

IVD

CE

Instrucciones de uso. 2009-11

EN	EXPLANATION OF SYMBOLS
BG	ОБЯСНЕНИЕ НА СИМВОЛИТЕ
CS	VÝZNAM SYMBOLŮ
DA	SYMBOLFORKLARING
DE	ERKLÄRUNG DER SYMBOLE
EL	ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES	SIGNIFICADO DE LOS SÍMBOLOS
ET	SÜMBOLITE SELGITUS
FR	EXPLICATION DES SYMBOLES
HR	OBJAŠNJENJE SIMBOLA
HU	JELMAGYARÁZAT
IT	SPIEGAZIONE DEI SIMBOLI
LT	SIMBOLIŲ PAAIŠKINIMAI
LV	SIMBOLU SKAIDROJUMS
NL	VERKLARING DER SYMBOLEN
NO	SYMBOLFORKLARING
PL	OBJAŚNIENIE SYMBOLI
PT	EXPLICAÇÃO DOS SÍMBOLOS
RO	SEMNIȚAȚIA SIMBOLURILOR
RU	ОБОЗНАЧЕНИЯ
SE	SYMBOLFÖRKLARING
SK	VÝZNAM SYMBOLOV
SL	RAZLAGA SIMBOLOV
SR	OBJAŠNJENJE SIMBOLA
TR	SEMBOLLERİN AÇIKLAMALARI

Use By/Годно до/Použitelné do/  
Holdbar til/Verwendbar bis/  
Ημερομηνία λήξης/Fecha  
de caducidad/Kölblük kuni/  
Utiliser jusque/Rok valjanosti/  
Felhasználható/Utilizzare entro/  
Sunautoti iki/Izlietot līdz/Houdbaar  
tot/Brukes innen/Użyć przed/  
Prazo de validade/Expirã la/  
Использовать до/Använd före/  
Použite né do/ Uporabno do/  
Upotrebljivo do/Son Kullanna Tarihi

LOT

Batch code/Номер на партида/  
Číslo šarže/Lotnummer/  
Chargenbezeichnung/Αριθμός  
Παρτίδας/Código de lote/Partii  
kood/Code du lot/Kod serije/  
Sarzsám/Codice del lotto/  
Partijos kodus/Partijas kods/Lot  
nummer/Partikode/Kod partii/  
Código do lote/Număr de lot/  
Номер лота/Lotnummer/Číslo  
šarže/Številka serije/Kod partije/  
Parti Kodu

Date of manufacture/Дата на производство/Datum výroby/  
Produktionsdato/Herstellungsdatum/  
Ημερομηνία παραγωγής/Fecha de fabricación/Valmistamise kuupäev/  
Date de fabrication/Datum proizvodnje/  
Gyártási idő/Data di produzione/  
Pagaminimo data/Ražošanas datums/  
Productiedatum/Fremstillingsdato/  
Data produkcji/Data de fabrico/Data fabricației/Дата производства/  
Tillverkningsdatum/Dátum výroby/Datum izdelave/Datum proizvodnje/Üretim tarihi

Temperature limitation/  
Температурни граници/  
Теплотни омеzeи/  
Temperaturbegränsning/  
Temperaturbegrenzung/  
Περιορισμοί θερμοκρασίας/  
Limites de temperatura/  
Temperatuuri piirang/  
Limite de température/  
Temperaturno ograničenje/  
Hőmérsékletre vonatkozó korlátozás/  
Limiti di temperatura/  
Temperatūriniai apribojimai/  
Temperatūras ierobežojums/  
Temperaturbepierking/  
Temperaturbegrensinger/  
Temperaturey graniczne/  
Limite de temperatura/  
Limite de temperatură/  
Температурный режим/  
Temperaturbegränsning/  
Teplotné obmedzenie  
Omejitve temperature/  
Temperaturno ograničenje/  
Sıcaklık sınırlaması/

## IVD

In Vitro Diagnostic Medical Device/  
Медицински уред за диагностика  
ин витро/Лéкаřský přístroj pro diagnostiku in vitro/Medicinsk udstyr til in vitro-diagnostik/In-vitro-Diagnostikum/  
Ιατροτεχνολογικό προϊόν για διάγνωση  
In Vitro/Dispositivo médico para diagnóstico in vitro/In vitro diagnostiline meditsiineasead/Dispositif médical de diagnostic in vitro/Diagnostički medicinski uređaj In Vitro/In vitro orvosdiagnosztikai eszköz/Dispositivo medico per test diagnostici in vitro/In Vitro Diagnostinė Medicinos Priemonė/  
Medicinska ierīce in vitro diagnostikai/  
In vitro-diagnostisch medisch instrument/  
In vitro diagnostisk medisinsk utstyr/  
Wyrób medyczny do diagnostyki in vitro/  
Dispositivo Médico de Diagnóstico In Vitro/Dispozitiv medical pentru diagnostic in vitro/Только для диагностики In Vitro/Endast för in vitro-diagnostik/  
Zdravotnička pomôcka na diagnostiku in vitro/In vitro diagnostični pripomoček/  
Diagnostički medicinski uređaj In Vitro/<96> testleri için yeterlilik içerir

Contains sufficient for <96> tests/Съдържа достатъчно количество за тестове <96>/Lze použít pro <96> testů/Ineholder tilstrækkeligt/Inhalt ausreichend für <96> Prüfungen/Περεχόμενο επαρκές για «96» εξετάσεις/Contenido suficiente para <96> ensayos/Kogusest piisab <96> testi läbiviimiseks/Contenu suffisant pour «96» tests/Sadržaj dovoljno za <96> testova/A doboz tartalma <96> vizsgálat elvégzéséhez elegendő/Contenuto sufficiente per «96» saggi/Turiny's skirtas atlikti <96> tyrimus/Saturs pietiekams <96> testiem/Inhoud voldoende voor «96» testen/til «96» test/ Tilstrækkelig innhold for <96> prøver/  
Wystarczy na wykonanie <96> testów/  
Conteúdo suficiente para «96» ensaios/  
Conținut suficient pentru 96 de teste/  
Содержит достаточные количества для «96» определений/Innehåller tillräckligt till «96» antal tester/Obsah postačuje na tento počet testov: <96>/Vsebina zadostuje za <96> testov/Sadržina dovoljna za <96> testova/<96> testleri için yeterlilik içerir

## REF

Catalogue number/Каталожен номер/  
Katalogové číslo/Katalognummer/  
Bestellnummer/Αριθμός καταλόγου/  
Número de catálogo/Katalogi number/  
Numéro de catalogue/Kataloški broj/  
Katalógusszám/Numero di catalogo/  
Katalogo numeris/Numurs katalogā/  
Catalogusnummer/Katalognummer/  
Numer katalogowy/Número do catálogo/  
Număr de catalog/Номер по каталогу/  
Produktnummer/Katalógové číslo/  
Kataloška številka/Kataloški broj/  
Katalog numarası

Consult Instructions for Use/  
Прочетете инструкцията за  
употреба/Konzultujte s návodem  
k použití/Se brugsanvisning/Siehe  
Gebrauchsanweisung/Συμβουλευτείτε  
της Οδηγίες σχετικά με τη χρήση/  
Consulte las instrucciones de uso/  
Vt kasutusjuhendit/Consulter le mode  
d'emploi/Pročítajte upute za uporabu/  
Olvassa el a használati utasítást/  
Consultare le istruzioni per l'uso/Dél  
naudojimo žiūrėkite instrukcijas/Izlasiet  
lietošanas instrukciju/Raadpleeg de  
instructies voor gebruik/Les instruksene  
for bruk/Sprawdzić w instrukcji użycia/  
Consulte as Instruções de Utilização/  
Consultați instrucțiunile de utilizare/  
Обратитесь к инструкции по  
применению/Se bruksanvisning/  
Prečítajte si návod na používanie/  
Pročítajte uputstvo za upotrebu/  
Kullanım Talimatlarını Bakınız

## CONT

Contents of kit/Съдържание на набора/  
Obsah sady/Kittets indhold/Inhalt des  
Kits/Περιεχόμενα του κιτ/Contenido  
del kit/Komplekt sisaldab/Contenu du  
kit/Sadržaj opreme/A készlet tartalma/  
Contenuto del kit/Rinkinio turinys/  
Komplekta saturs/Inhoud van de set/  
Settets innhold/Zawartość zestawu/  
Conteúdo do kit/Conținutul setului/  
Компоненты набора/Kit innehåll/  
Obsah súpravy/Vsebina kompleta/Sadržaj  
opreme/Kitin içindekiler

Biological risks/Биологическа  
опасност/Biológická rizika/Biologisk  
fare/Biologische Gefahren/Biológikoi  
kínðuvon/Riesgos biológicos/  
Biolooigilised ohud/Risques biologiques/  
Biolóškli rizici/Biológiai kockázatok/Rischi  
biologici/Biologinis pavojus/Biológiskais  
risks/Biologische risico's/Biologiske  
risikoer/Zagroženie biologiczne/Riscos  
biológicos/ Biologisk risk/Pericole  
biologice/Биологическая опасность/  
Biologický rizikové/Biológické riziká/  
Biolóškli rizici/Biyolojik riskler

## ORIG HUM

Human/C човешки производ/Lidské/  
Human/Human/δείγματα αναφοράς/  
Humano/Inimpãritolu/Humaine/Ljudskog  
porjekla/Humãn/Origine Umana/  
Žmogaus kilmės/Cilvēku izcelsmes/  
Human/Menneske/Ludzka/Humano/  
Origine umãnã/Человеческого  
происхождения/Human/Ludské/  
Humanega izvora/Ljudskog porekla/İnsan

## ORIG MOU

From mouse/C миши производ/Myši/  
Fra mus/der Maus/από ποντίκι/de raton/  
Hiirtelt/De souris/Mišijeg porjekla/  
Egérbdli/Murino/Pelès kilmės/No peles/  
Van muizen/Fra mus/Mysia/Do rato/De  
la șoareci/Мышиного происхождения/  
Från mus/Myšije/Mišjega izvora/Mišijeg  
porekla/Fareden

## ORIG BOV

Bovine/C говежди производ/  
Hovēži/Bovin/Rind/από βοοειδή/  
Bovino/Veistelt/Bovine/Rogate stoke/  
Szarvasmarha/Bovina/Jaučio/No  
liellopa/Bovien/Bovina/Wolowy/Bovino/  
Origine bovinã/крупного рогатого  
скота/Från ko/Hovädzie/Rogaveja  
izvora/Rogate krupne stoke/Bovin

Reconstitute with/Разтваряне с/  
Rozfeđe pomoci/Rekonstitutes med/  
Rekonstituieren mit/Ανασύσταση με/  
Reconstituir con/Lahjendamine/  
Reconstituer avec/Rekonstituiraite s/  
Feloldáshoz/Ricostituire con/LT/Atškaidīt  
ar/Reconstituie me/Rekonstituerees  
med/Odtworzyć za pomocą/Reconstituir  
com/A se reconstitui cu/Растворить в/  
Rekonstituera med/Rozriedte pomocou/  
Rekonstituiraite z/s/ Ponadto formiranje  
sa/Yeniden oluřturulur

Manufacturer/Производитель/Výrobce/  
Producent/Hersteller/Κτασκευαστής/  
Fabricante/Tootja/Fabricant/Proizvođač/  
Gyártó/Fabbricante/Gamintojas/  
Ražotājs/Fabrikant/Produsent/  
Producent/Fabricante/Producător/  
Производитель/Тилverkare/ Výrobca/  
Izdevalavec/Proizvođač/Üretici

# CanAg NSE EIA

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1. Paus E. and Nustad K., (1989) Immunoradiometric Assay for  $\alpha\gamma$ - and  $\gamma\gamma$ -Enolase (Neuron-Specific Enolase), with Use of Monoclonal Antibodies and Magnetizable Polymer Particles. *Clin. Chem.* 35: 2034-2038.
2. Dahlén U., Karlsson B., Nilsson O. and Uhl W., (1995) Development of an Enzyme Immunoassay, NSE-Enzymun Test For Determination of Neuron-Specific Enolase. XXIII International Society for Oncodevelopmental Biology and Medicine, Montréal, Québec .
3. Cooper E.H., (1994) Neuron-specific enolase. *The International Journal of Biological Markers* 9(4):205-10.
4. Cooper E.H., Pritchard J., Bailey C.C. and Ninane J., (1987) Serum neuron-specific enolase in children's cancer. *Br. J. Cancer* 56: 65-67.
5. Schneider, P. M. et al., (2002) Lung Cancer. In "Tumor markers, Physiology, Pathobiology, Technology and Clinical Applications" Eds. Diamandis E. P. et al., AACCC Press, Washington pp 287-303.
6. Bonner J. A., Sloan JA., Rowland KM., Klee GG., Kugler JW., Mailliard JA., Wiesenfeld M., Krook JE., Maksymiuk AW., Shaw EG., Marks RS and Perez EA., (2000) Significance of Neuron-specific Enolase Levels before and during Therapy for Small Cell Lung Cancer. *Clinical Cancer Research* 6: 597-601.
7. Pählman S., Esscher T., Bergvall P. and Odelstad L., (1984) Purification and characterization of human neuron-specific enolase: Radioimmunoassay development. *Tumor Biol.* 5: 127-139.
8. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
9. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).

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**CanAg®** es una marca registrada de Fujirebio Diagnostics AB
**Fujirebio Diagnostics AB**